



# SPAC1420.01c Antibody

|                     |                                                                                                      |
|---------------------|------------------------------------------------------------------------------------------------------|
| Product Code        | CSB-PA892050XA01SXV                                                                                  |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                        |
| Uniprot No.         | Q9UTM5                                                                                               |
| Immunogen           | Recombinant Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast) SPAC1420.01c protein |
| Raised In           | Rabbit                                                                                               |
| Species Reactivity  | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)                                  |
| Tested Applications | ELISA, WB (ensure identification of antigen)                                                         |
| Form                | Liquid                                                                                               |
| Conjugate           | Non-conjugated                                                                                       |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                     |
| Purification Method | Antigen Affinity Purified                                                                            |
| Isotype             | IgG                                                                                                  |
| Clonality           | Polyclonal                                                                                           |
| Product Type        | Polyclonal Antibody                                                                                  |
| Immunogen Species   | Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)                                  |
| Lead Time           | Made-to-order (14-16 weeks)                                                                          |
| Target Names        | SPAC1420.01c                                                                                         |
| Usage               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                          |