



NADH dehydrogenase [ubiquinone] iron-sulfur protein 8, mitochondrial Antibody

Product Code	CSB-PA305275XA01FIG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80269
Immunogen	Recombinant Solanum tuberosum (Potato) NADH dehydrogenase [ubiquinone] iron-sulfur protein 8, mitochondrial protein
Raised In	Rabbit
Species Reactivity	Solanum tuberosum (Potato)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Solanum tuberosum (Potato)
Lead Time	Made-to-order (14-16 weeks)