



StPat14K07.03 Antibody

Product Code	CSB-PA643598XA01FIG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2MY60
Immunogen	Recombinant Solanum tuberosum (Potato) StPat14K07.03 protein
Raised In	Rabbit
Species Reactivity	Solanum tuberosum (Potato)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Solanum tuberosum (Potato)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	StPat14K07.03