



Pyruvate dehydrogenase E1 component subunit beta, mitochondrial Antibody

Product Code	CSB-PA305505XA01FIG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P81419
Immunogen	Recombinant Solanum tuberosum (Potato) Pyruvate dehydrogenase E1 component subunit beta, mitochondrial protein
Raised In	Rabbit
Species Reactivity	Solanum tuberosum (Potato)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Solanum tuberosum (Potato)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.