



PETN Antibody

Product Code	CSB-PA652789XA01FIG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2VEI3
Immunogen	Recombinant Solanum tuberosum (Potato) PETN protein
Raised In	Rabbit
Species Reactivity	Solanum tuberosum (Potato)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Species	Solanum tuberosum (Potato)
Gene Names	petN