



QCR9 Antibody

Product Code	CSB-PA343185XA01FIG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46270
Immunogen	Recombinant Solanum tuberosum (Potato) QCR9 protein
Raised In	Rabbit
Species Reactivity	Solanum tuberosum (Potato)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Species	Solanum tuberosum (Potato)