



rpoC2 Antibody

Product Code	CSB-PA647294XA01FIG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2VEI6
Immunogen	Recombinant Solanum tuberosum (Potato) rpoC2 protein
Raised In	Rabbit
Species Reactivity	Solanum tuberosum (Potato)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Solanum tuberosum (Potato)
Lead Time	Made-to-order (14-16 weeks)
Target Names	rpoC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.