



Long neurotoxin MS2 Antibody

Product Code	CSB-PA308464XA01MGQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P86096
Immunogen	Recombinant Micrurus surinamensis (Surinam coral snake) Long neurotoxin MS2 protein
Raised In	Rabbit
Species Reactivity	Micrurus surinamensis (Surinam coral snake)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Micrurus surinamensis (Surinam coral snake)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.