



afp18 Antibody

Product Code	CSB-PA270774XA01GIA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C4ULG3
Immunogen	Recombinant Yersinia ruckeri serotype O1 (strain ATCC 29473 / DSM 18506 / JCM 15110 / NCIMB 2194 / NCTC 12986 / 2396-61) afp18 protein
Raised In	Rabbit
Species Reactivity	Yersinia ruckeri serotype O1 (strain ATCC 29473 / DSM 18506 / JCM 15110 / NCIMB 2194 / NCTC 12986 / 2396-61)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Yersinia ruckeri serotype O1 (strain ATCC 29473 / DSM 18506 / JCM 15110 / NCIMB 2194 / NCTC 12986 / 2396-61)
Lead Time	Made-to-order (14-16 weeks)
Target Names	afp18
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.