



conoCAP Antibody

Product Code	CSB-PA517182XA01DWL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	E3PQQ8
Immunogen	Recombinant Conus villepinii (Villepin's cone) conoCAP protein
Raised In	Rabbit
Species Reactivity	Conus villepinii (Villepin's cone)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Conus villepinii (Villepin's cone)
Lead Time	Made-to-order (14-16 weeks)
Target Names	conoCAP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.