



cas9 Antibody

Product Code	CSB-PA166726XA01SMT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99ZW2
Immunogen	Recombinant Streptococcus pyogenes serotype M1 cas9 protein
Raised In	Rabbit
Species Reactivity	Streptococcus pyogenes serotype M1
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Streptococcus pyogenes serotype M1
Lead Time	Made-to-order (14-16 weeks)
Gene Names	cas9