



CK054 Antibody

Product Code	CSB-PA764679XA01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NWE0
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) CK054 protein
Raised In	Rabbit
Species Reactivity	Danio rerio (Zebrafish) (Brachydanio rerio)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugate
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Gene Names	zgc:85789