



CNPY1 Antibody

Product Code	CSB-PA646601XA01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2L6L1
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) CNPY1 protein
Raised In	Rabbit
Species Reactivity	Danio rerio (Zebrafish) (Brachydanio rerio)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Gene Names	cnpy1