



UPF0729 protein C18orf32 homolog Antibody

Product Code	CSB-PA252631XA01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C8Y2
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) UPF0729 protein C18orf32 homolog protein
Raised In	Rabbit
Species Reactivity	Danio rerio (Zebrafish) (Brachydanio rerio)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Lead Time	Made-to-order (14-16 weeks)