



UPF0545 protein C22orf39 homolog Antibody

Product Code	CSB-PA609514XA01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q0P3X7
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) UPF0545 protein C22orf39 homolog protein
Raised In	Rabbit
Species Reactivity	Danio rerio (Zebrafish) (Brachydanio rerio)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.