



DTWD1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA736465XA01DIL |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q5XJ56 |
| Immunogen | Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) DTWD1 protein |
| Raised In | Rabbit |
| Species Reactivity | Danio rerio (Zebrafish) (Brachydanio rerio) |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugate |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Protein A/G |
| Isotype | IgG |
| Clonality | Polyclonal |
| Product Type | Polyclonal Antibody |
| Species | Danio rerio (Zebrafish) (Brachydanio rerio) |
| Gene Names | dtwd1 |