



EPA4A Antibody

Product Code	CSB-PA835524XA01XBE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q91845
Immunogen	Recombinant Xenopus laevis (African clawed frog) EPA4A protein
Raised In	Rabbit
Species Reactivity	Xenopus laevis (African clawed frog)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Species	Xenopus laevis (African clawed frog)
Gene Names	epha4-a