



foxg1 Antibody

Product Code	CSB-PA897423XA01XBE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9YHC5
Immunogen	Recombinant Xenopus laevis (African clawed frog) foxg1 protein
Raised In	Rabbit
Species Reactivity	Xenopus laevis (African clawed frog)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Xenopus laevis (African clawed frog)
Lead Time	Made-to-order (14-16 weeks)
Target Names	foxg1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.