



fen1 Antibody

Product Code	CSB-PA686312XA01XAS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5I4H3
Immunogen	Recombinant Xiphophorus maculatus (Southern platyfish) (Platypoecilus maculatus) fen1 protein
Raised In	Rabbit
Species Reactivity	Xiphophorus maculatus (Southern platyfish) (Platypoecilus maculatus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Xiphophorus maculatus (Southern platyfish) (Platypoecilus maculatus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	fen1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.