



foxi1e Antibody

Product Code	CSB-PA773654XA01XBE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7ZYQ0
Immunogen	Recombinant Xenopus laevis (African clawed frog) foxi1e protein
Raised In	Rabbit
Species Reactivity	Xenopus laevis (African clawed frog)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Xenopus laevis (African clawed frog)
Lead Time	Made-to-order (14-16 weeks)
Target Names	foxi1e
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.