



gnb2l1 Antibody

Product Code	CSB-PA009605XA01OET
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O42249
Immunogen	Recombinant Oreochromis niloticus (Nile tilapia) (Tilapia nilotica) gnb2l1 protein
Raised In	Rabbit
Species Reactivity	Oreochromis niloticus (Nile tilapia) (Tilapia nilotica)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Oreochromis niloticus (Nile tilapia) (Tilapia nilotica)
Lead Time	Made-to-order (14-16 weeks)
Target Names	gnb2l1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.