



gpr183b Antibody

Product Code	CSB-PA535153XA01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B0UXR0
Immunogen	Synthetic Peptide (1-25aa)
Species Reactivity	Danio rerio (Zebrafish) (Brachydanio rerio)
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Target Names	gpr183b
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.