



ihh Antibody

Product Code	CSB-PA011569XA01DJX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O13243
Immunogen	Recombinant Devario devario (Bengal danio) (Danio devario) ihh protein
Raised In	Rabbit
Species Reactivity	Devario devario (Bengal danio) (Danio devario)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Devario devario (Bengal danio) (Danio devario)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ihh
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.