



ihh Antibody

Product Code	CSB-PA011569XA01TIJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P79866
Immunogen	Recombinant Trigonostigma heteromorpha (Harlequin rasbora) (Rasbora heteromorpha) ihh protein
Raised In	Rabbit
Species Reactivity	Trigonostigma heteromorpha (Harlequin rasbora) (Rasbora heteromorpha)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Trigonostigma heteromorpha (Harlequin rasbora) (Rasbora heteromorpha)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ihh
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.