



Con-Ins T1A Antibody

Product Code	CSB-PA553007XA01CPL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0A0B5ADU4
Immunogen	Recombinant Conus tulipa (Fish-hunting cone snail) (Tulip cone) Con-Ins T1A protein
Raised In	Rabbit
Species Reactivity	Conus tulipa (Fish-hunting cone snail) (Tulip cone)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Conus tulipa (Fish-hunting cone snail) (Tulip cone)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.