



Ihx5 Antibody

Product Code	CSB-PA012923XA01XBE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37137
Immunogen	Recombinant Xenopus laevis (African clawed frog) Ihx5 protein
Raised In	Rabbit
Species Reactivity	Xenopus laevis (African clawed frog)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Xenopus laevis (African clawed frog)
Lead Time	Made-to-order (14-16 weeks)
Target Names	Ihx5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.