



ngb Antibody

Product Code	CSB-PA015776XA01TCA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q90W04
Immunogen	Recombinant Tetraodon nigroviridis (Spotted green pufferfish) (Chelonodon nigroviridis) ngb protein
Raised In	Rabbit
Species Reactivity	Tetraodon nigroviridis (Spotted green pufferfish) (Chelonodon nigroviridis)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Tetraodon nigroviridis (Spotted green pufferfish) (Chelonodon nigroviridis)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	ngb