



bglapBGP Antibody

Product Code	CSB-PA556032XA01STX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q1EG28
Immunogen	Recombinant Solea senegalensis (Senegalese sole) bglapBGP protein
Raised In	Rabbit
Species Reactivity	Solea senegalensis (Senegalese sole)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Solea senegalensis (Senegalese sole)
Lead Time	Made-to-order (14-16 weeks)
Target Names	bglapBGP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.