



pkd2 Antibody

Product Code	CSB-PA439269XA01GID
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	H2LRU7
Immunogen	Recombinant Oryzias latipes (Japanese rice fish) (Japanese killifish) pkd2 protein
Raised In	Rabbit
Species Reactivity	Oryzias latipes (Japanese rice fish) (Japanese killifish)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Oryzias latipes (Japanese rice fish) (Japanese killifish)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	pkd2