



slc39a1 Antibody

Product Code	CSB-PA764841XA01TKY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6QQT1
Immunogen	Recombinant Takifugu rubripes (Japanese pufferfish) (Fugu rubripes) slc39a1 protein
Raised In	Rabbit
Species Reactivity	Takifugu rubripes (Japanese pufferfish) (Fugu rubripes)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Takifugu rubripes (Japanese pufferfish) (Fugu rubripes)
Lead Time	Made-to-order (14-16 weeks)
Target Names	slc39a1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.