



shh Antibody

Product Code	CSB-PA021266XA01DIJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O13235
Immunogen	Recombinant Danio aff. albolineatus shh protein
Raised In	Rabbit
Species Reactivity	Danio aff. albolineatus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Danio aff. albolineatus
Lead Time	Made-to-order (14-16 weeks)
Gene Names	shh