



shh Antibody

Product Code	CSB-PA021266XA01DIK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O13238
Immunogen	Recombinant Danio aff. tweediei shh protein
Raised In	Rabbit
Species Reactivity	Danio aff. tweediei
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Danio aff. tweediei
Lead Time	Made-to-order (14-16 weeks)
Gene Names	shh