



shh Antibody

Product Code	CSB-PA021266XA01GIB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P79709
Immunogen	Recombinant Danio kerri (Blue danio) (Brachydanio kerri) shh protein
Raised In	Rabbit
Species Reactivity	Danio kerri (Blue danio) (Brachydanio kerri)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Danio kerri (Blue danio) (Brachydanio kerri)
Lead Time	Made-to-order (14-16 weeks)
Target Names	shh
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.