



shh Antibody

Product Code	CSB-PA021266XA01GIC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P79850
Immunogen	Recombinant Puntigrus tetrazona (Sumatra barb) (Puntius tetrazona) shh protein
Raised In	Rabbit
Species Reactivity	Puntigrus tetrazona (Sumatra barb) (Puntius tetrazona)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Puntigrus tetrazona (Sumatra barb) (Puntius tetrazona)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	shh