



shh Antibody

Product Code	CSB-PA021266XA01GIE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P79869
Immunogen	Recombinant Rasbora paviana (Sidestripe rasbora) shh protein
Raised In	Rabbit
Species Reactivity	Rasbora paviana (Sidestripe rasbora)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rasbora paviana (Sidestripe rasbora)
Lead Time	Made-to-order (14-16 weeks)
Target Names	shh
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.