



shh Antibody

Product Code	CSB-PA021266XA01TAR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P79915
Immunogen	Recombinant Tanichthys albonubes (White cloud mountain minnow) shh protein
Raised In	Rabbit
Species Reactivity	Tanichthys albonubes (White cloud mountain minnow)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Tanichthys albonubes (White cloud mountain minnow)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	shh