



shh Antibody

Product Code	CSB-PA021266XA01DJY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O13247
Immunogen	Recombinant Devario malabaricus (Malabar danio) (Danio malabaricus) shh protein
Raised In	Rabbit
Species Reactivity	Devario malabaricus (Malabar danio) (Danio malabaricus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Devario malabaricus (Malabar danio) (Danio malabaricus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	shh
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.