



snap25a Antibody

Product Code	CSB-PA334206XA01DQG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36977
Immunogen	Recombinant Carassius auratus (Goldfish) snap25a protein
Raised In	Rabbit
Species Reactivity	Carassius auratus (Goldfish)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Carassius auratus (Goldfish)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	snap25a