



# syndig1l Antibody

|                     |                                                                                  |
|---------------------|----------------------------------------------------------------------------------|
| Product Code        | CSB-PA729231XA01DIL                                                              |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| Uniprot No.         | Q5TZE2                                                                           |
| Immunogen           | Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) syndig1l protein         |
| Raised In           | Rabbit                                                                           |
| Species Reactivity  | Danio rerio (Zebrafish) (Brachydanio rerio)                                      |
| Tested Applications | ELISA, WB (ensure identification of antigen)                                     |
| Form                | Liquid                                                                           |
| Conjugate           | Non-conjugated                                                                   |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified                                                        |
| Isotype             | IgG                                                                              |
| Clonality           | Polyclonal                                                                       |
| Product Type        | Polyclonal Antibody                                                              |
| Immunogen Species   | Danio rerio (Zebrafish) (Brachydanio rerio)                                      |
| Lead Time           | Made-to-order (14-16 weeks)                                                      |
| Target Names        | syndig1l                                                                         |
| Usage               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.      |