



tbpl2 Antibody

Product Code	CSB-PA732985XA01XBE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5UE94
Immunogen	Recombinant Xenopus laevis (African clawed frog) tbpl2 protein
Raised In	Rabbit
Species Reactivity	Xenopus laevis (African clawed frog)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Xenopus laevis (African clawed frog)
Lead Time	Made-to-order (14-16 weeks)
Target Names	tbpl2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.