



TRADD Antibody

Product Code	CSB-PA867432XA01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9I9N5
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) TRADD protein
Raised In	Rabbit
Species Reactivity	Danio rerio (Zebrafish) (Brachydanio rerio)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Lead Time	Made-to-order
Gene Names	tradd