



tnfaip8l2a Antibody

Product Code	CSB-PA624751XA01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q1ECV8
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) tnfaip8l2a protein
Raised In	Rabbit
Species Reactivity	Danio rerio (Zebrafish) (Brachydanio rerio)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Lead Time	Made-to-order (14-16 weeks)
Target Names	tnfaip8l2a
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.