



twhh Antibody

Product Code	CSB-PA301286XA01GIC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P79855
Immunogen	Recombinant Puntigrus tetrazona (Sumatra barb) (Puntius tetrazona) twhh protein
Raised In	Rabbit
Species Reactivity	Puntigrus tetrazona (Sumatra barb) (Puntius tetrazona)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Puntigrus tetrazona (Sumatra barb) (Puntius tetrazona)
Lead Time	Made-to-order (14-16 weeks)
Target Names	twhh
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.