



U2AF1L4 Antibody

Product Code	CSB-PA667525XA01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3T127
Immunogen	Recombinant Bos taurus (Bovine) U2AF1L4 protein
Raised In	Rabbit
Species Reactivity	Bos taurus (Bovine)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bos taurus (Bovine)
Lead Time	Made-to-order (14-16 weeks)
Target Names	U2AF1L4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.