



U2af1I4 Antibody

Product Code	CSB-PA805251XA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8BGJ9
Immunogen	Recombinant Mus musculus (Mouse) U2af1I4 protein
Raised In	Rabbit
Species Reactivity	Mus musculus (Mouse)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Mus musculus (Mouse)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	U2af1I4