



twhh Antibody

Product Code	CSB-PA517523XA01DJZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O13253
Immunogen	Recombinant Devario pathirana (Barred danio) (Danio pathirana) twhh protein
Raised In	Rabbit
Species Reactivity	Devario pathirana (Barred danio) (Danio pathirana)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Devario pathirana (Barred danio) (Danio pathirana)
Lead Time	Made-to-order (14-16 weeks)
Target Names	twhh
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.