



twhh Antibody

Product Code	CSB-PA519261XA01FGX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O13226
Immunogen	Recombinant Pseudorasbora parva twhh protein
Raised In	Rabbit
Species Reactivity	Pseudorasbora parva
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pseudorasbora parva
Lead Time	Made-to-order (14-16 weeks)
Target Names	twhh
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.