



atp6v1c1a Antibody

Product Code	CSB-PA763084XA01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7T385
Immunogen	Recombinant Danio rerio (Zebrafish) (Brachydanio rerio) atp6v1c1a protein
Raised In	Rabbit
Species Reactivity	Danio rerio (Zebrafish) (Brachydanio rerio)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Lead Time	Made-to-order (14-16 weeks)
Target Names	atp6v1c1a
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.