





Recombinant Human N-terminal Xaa-Pro-Lys N-methyltransferase 1 (NTMT1) (Active)

| | |
|----------------------------|---|
| Product Code | CSB-EP880111HU1 |
| Abbreviation | Recombinant Human NTMT1 protein (Active) |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9BV86 |
| Form | Liquid or Lyophilized powder |
| Storage Buffer | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose. |
| Product Type | Recombinant Human N-terminal Xaa-Pro-Lys N-methyltransferase 1(NTMT1) |
| Immunogen Species | Homo sapiens (Human) |
| Biological Activity | Activity Measured by its binding ability in a functional ELISA. Immobilized CCL5 at 2 µg/ml can bind human NTMT1, the EC ₅₀ of human NTMT1 protein is 3.029-5.404 µg/ml. |
| Purity | ≥ 85% as determined by SDS-PAGE. |
| Sequence | TSEVIEDEKQFYKAKTYWKQIPPTVDGMLGGYGHISSIDINSSRKFLQRFLRE GPNKTGTSCALDCGAGIGRITKRLLLPLFREVDMDITEDFLVQAKTYLGEEGK RVRNYFCCGLQDFTPEPDSYDVIWIQWVIGHLTDQHAEFLRRCKGSLRPNGLI VIKDNMAQEGVILDDVDSSVCRDLDDVRRRIICSAGLSLLAEERQENLPPDEIYHVY SFALR |
| Research Area | Cell Biology |
| Source | E.coli |
| Target Names | NTMT1 |
| Expression Region | 2-223aa |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | N-terminal GST-tagged |
| Mol. Weight | 52.7 kDa |
| Protein Length | Full Length of Mature Protein |
| Image | |



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Activity

Measured by its binding ability in a functional ELISA. Immobilized CCL5 at 2 µg/ml can bind human NTMT1, the EC₅₀ of human NTMT1 protein is 3.029-5.404 µg/ml.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.