




# Recombinant Clostridioides difficile NAD-specific glutamate dehydrogenase (gluD)

<b>Product Code</b>	CSB-EP335284DUNa0
<b>Abbreviation</b>	Recombinant Clostridioides difficile gluD protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P27346
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Clostridioides difficile (Peptoclostridium difficile)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	MSGKDVNVFEMAQSQVKNACDKLGMEPAVYELLKEPMRVIEVSIPVKMDDGS IKTFKGFERSQHNDVAVGPTKGGIRFHQNVSRDEVKALSIWMTFKCSVTGIPYGG GKGGIIVDPSTLSQGELERLSRGYIDGIYKLIGEKVDVPAPDVNTNGQIMSWMV DEYNKLTGQSSIGVITGKPVVEFGSLGRATAATGFGVAVTAREAAKLGIDMKK AKIAVQQGIGNVGSYTVLNCEKLGGTVVAMAEWCKSEGSYAIYNENGLDGGQAM LDYMKEHGNLLNFPGAKRISLEEFWASDVIDIVIPAALENSITKEVAESIKAKLVC EAANGPTTPEADEVFAERGIVLTPDILTANAGGVTVSIFYFEWVQNLGYGYWSEEE VEQKEEIAMVKAFESIWKIKEEYNVTMREAAVMHSIKKVAEAMKLRGWY
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	gluD
<b>Expression Region</b>	1-421aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	50.0 kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.